

ATLAS™ FFN design features

Syncera
Powered by  smith&nephew

- 5° lateral offset for minimally invasive trochanteric entry
- 12° of built-in femoral neck ante-version
- Neck shaft angle of 130 degrees
- 2 M-L view distal locking options including one dynamic slot
- Nail diameters – 9, 10, 11.5 & 13mm
- Nail length range: - 30cm to 50cm
- Screws
 - 6.4 mm Recon screws (Length range - 65 to 125mm)
 - 5 mm screws (gold – from 20mm to 80mm in 2.5mm increments & 80mm to 110mm in 5mm increments) for 10, 11.5 & 13mm nails
 - 4.5 mm (green – from 25mm to 65mm in 5mm increments) for 9 mm nails
- End Cap: 0 mm, 5 mm & 10 mm



Not for distribution. For Limited Market Release only.